VECTEC

Manufactured by :

VECTEC


Bioparc - 03270 HAUTERIVE - FRANCE

Tél. : + 33 (0)4 70 58 41 70

Fax : + 33 (0)4 70 58 41 79

vectec@vectec.fr

REF AW16280

F	Planche A :	PROCÉDURE DE RETRAIT DU T'LIFT®	
GB	Plate A :	T'LIFT® WITHDRAWAL PROCEDURE	
D	Bildtafel A :	RÜCKZUGSVERFAHREN DES T'LIFT®	
E	Plancha A :	PROCEDIMIENTO DE RETIRADA DEL T'LIFT®	
I	Tavola A :	PROCEDURA DI RITIRO DEL T'LIFT®	
P	Quadro A :	PROCESSO DE RETRAÇÃO DO T'LIFT®	
Dk	Diagram A :	UDTAGNINGSPROCEDURE FOR T'LIFT®	
N	Plansje A :	UTTREKKINGSPROSEDYRE FOR T'LIFT®	
S	Bild A :	PROCEDUR FÖR ATT AVLAGSNA T'LIFT®	
		①	②
F	Coupe	Prise	Retrait
GB	Cut	Hold	Withdrawal
D	Schnitt	Griff	Rückzug
E	Corte	Sujeción	Retirada
I	Taglio	Presa	Ritiro
P	Corte	Tomada	Retração
Dk	Skæring	Greb	Udtagning
N	Kutt	Grip	Trekk ut
S	Klipp	Fatta tag	Avlägsna

F

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CAUTION

Please carefully read the present Instructions prior to using the VECTEC single-use T'Lift® laparoscopic retraction system. Failure to comply with the following recommendations is liable to entail serious consequences. Endoscopic operations should only be performed by qualified doctors having been specifically trained for the purpose and having gained sufficient experience of use. All surgeons performing endoscopic operations should be aware of the possible risks and complications involved using said surgical technique (notably risks for the patient and in manipulating instruments connected to an external power source, such as electrosurgery, ultrasound, laser, gas and liquid pressure). Pursuant to legislation in force, all patients having a right to information, adequate details of the operation should be given.

WARNINGS

- ⚠ **These instructions contain information on the VECTEC single-use T'Lift® only, as well as on how to use it. On no account should they be considered as an instructions' manual for using laparoscopic retraction systems in general.**
- ⚠ **The VECTEC T'Lift® is designed for single use only. On no account may it be used a second time over, whatever the reason or the means.**
- ⚠ **The device is delivered in sterile packaging. The T'Lift® should never be resterilised, regardless of whether it has been used or not.**
- ⚠ **Any reuse of this instrument shall affect its original sterility, may induce abnormal functioning or may lead to a cross-contamination in patients. The manufacturer may not be held liable in such a case. Reuse or resterilisation may also compromise the integrity of the instrument and use thereof, possibly causing injury, sickness or even death of the patient.**
- ⚠ **Warnings specific to a particular manipulation are stated herein, throughout the present Instructions (in bold type and preceded by ⚠)**
- ⚠ **Responsibility with respect to the storage, opening, appropriate use, handling, insertion and destruction of the VECTEC T'Lift® lies solely with the operator. On no account may VECTEC be held liable in the case of failure by the operator to observe proper surgical practices for the storage, the opening, the appropriate use, the manipulation, the insertion and the destruction of the VECTEC T'Lift®.**

INDICATIONS

The VECTEC T'Lift® is an instrument exclusively designed for endoscopic surgery. It should only be used by qualified practitioners having gained due experience. The T'Lift® is designed to enable the suspension and the presentation of anatomical parts throughout endoscopic surgery.

CONTRAINDICATIONS

Use of the VECTEC T'Lift® is strictly forbidden whenever endoscopic surgery is not advisable in a patient. On no account may VECTEC be held liable for assessment of the contraindication in relation to endoscopic surgery.

DESCRIPTION OF THE INSTRUMENT

The laparoscopic retraction system is packed along with:

- 2 T'Lifts®,
- 2 locking clamps,
- 1 applicator (body and needle).

INSTRUCTIONS FOR USE

- 1- Make sure that the packaging (blister) guaranteeing sterile conditions for the medical instrument is neither open nor damaged.
- 2- Check the use-by date printed on the packaging.
- 3- Open the packaging according to the proper practices for zones of operation and place the 5 elements on the operating zone.

⚠ **The VECTEC laparoscopic retraction system and its packaging should be destroyed as of when sterility may no longer be guaranteed, notably in the following cases: expiry of the use-by date on the packaging; packaging severely crushed or clearly damaged; damp marks visible on the packaging; perforation of the blister (shell or protective cover).**

- 4- Check:
 - the good state of the laparoscopic retraction system;
 - the existence of all components (T'Lift®, clamps and applicator).

⚠ **A damaged laparoscopic retraction system should never be used!**

5- Before starting surgery, make sure that all endoscopic instruments may be fitted to the VECTEC T'Lift®. Before inserting the VECTEC T'Lift®, the operating area, i.e. the abdominal wall, should be prepared as usual using proper surgical practices.

⚠ **On no account may VECTEC be held liable in the case of failure by the operator to observe proper surgical practices when verifying whether the endoscopic instruments are compatible with VECTEC T'Lift®.**

6- Remove the protection tube from the applicator; the user should pay strict attention not to injure himself when handling the applicator (sharp object).

7- Insert the distal part of the T'Lift® into the main body of the applicator and gently push until reaching the first marker.

8- Grasp the T'Lift® fitted with the applicator and perforate the abdominal wall, using proper surgical practices on the one hand and checking the trajectory and exit point of the needle in the abdominal cavity under endoscopy on the other.

9- Under endoscopy, grasp the anatomical part to be suspended, position it on the needle of the applicator and guide it through the anatomical part.

10- Hold the body of the applicator and gently push the T'Lift® shaft until reaching the second marker.

11- Under endoscopy, check the correct opening of the T'Lift® beneath the anatomical part chosen to be suspended.

12- Do not pull on the T'Lift® before having removed the applicator in order to prevent damaging it.

13- Hold the T'Lift® in position using an atraumatic clamp and remove the applicator from the outside, without releasing the anatomical part.

14- Place a clamp on the T'Lift® and block it in the desired position in order to correctly suspend the anatomical part. Do not "pull" the part solely with the T'Lift®, but guide it using an atraumatic clamp.

⚠ **Be sure to follow the regular surgical techniques according to proper surgical practices in order to protect the vessels and organs from any possible injury.**

15- Do not exert any undue tension on the anatomical part and check the quality of the tissues suspended.

⚠ **On no account may VECTEC be held responsible in the case of failure by the operator to comply with the recommendations for surgical techniques.**

⚠ **Do not use any electrosurgical instruments directly in contact with the T'Lift®: risk of damage to the T'Lift®.**

16- When suspension is no longer required, follow the T'Lift® withdrawal procedure (see plate A): cut the T'Lift® on the outside of the abdominal cavity (see figure ①) beneath the clamp, thereby enabling to easily identify the part of the T'Lift® inside the abdominal cavity. Then extract this part using a trocar. The T'Lift® should be held in place by its distal end (the largest and most curved part of the T'Lift®: see figure ②) using a clamp. Gently pull out the T'Lift® (see figure ③) under endoscopic control, in order not to damage the tissues in the abdominal cavity or to sever any part of the T'Lift®.

17- Check immediately the two parts of the T'Lift® on the operation zone and make sure that the instrument is complete.

18- Repeat the above actions (from point 7 to point 17) for the 2nd T'Lift® and the 2nd clamp.

Upon completion of the procedure, the T'Lift® should be disposed of as operating-theatre waste:

- All T'Lifts®, clamps and protective tubes should be placed in special containers reserved for soiled instruments,
- The applicator should be placed in the container reserved for sharp, traumatizing or perforating objects.

WARNING

Never use and thus destroy the sterile packaging and its contents in the following cases:

- Packaging severely crushed or clearly damaged.
- Perforation of the packaging (protective cover and/or hard shell).

STORAGE CONDITIONS

- Storage temperature: avoid any prolonged exposure to low and/or high temperatures.
- Store in an environment adapted to sterile medical devices.
- For transportation, all the above points should be applied.

Information for the United States (USA):

- US Federal Law restricts the sale of this instrument to that by doctors or upon prescription only.
- Sterilised with ethylene oxide.
- Single use only - Not to be reused.

In the event of any incident using this product (or medical device vigilance), please contact VECTEC at the following address: regulatory@vectec.fr

Initial CE marking: DECEMBER 2010.

Last update of these instructions: NOVEMBER 2012.

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