



## DISPOSABLE INSUFFLATION NEEDLES

Reference	Length	Diameter
AIG2 120	120 mm	2,1 mm
AIG2 150	150 mm	2,1 mm

The ergonomic body provides a comfortable hold and a clear feeling during the introduction.

The atraumatic blunt tip avoids the risk of perforation of intra and retroperitoneal organs.

The movement of a red-colored indicator inside the transparent body indicates the passing through of the various layers of the abdominal wall.

The specific sharpening of the needle guarantees a perfect perforation.

Single use permanently provides an optimal level of quality.

The needle is fitted with a stopcock and an integrated luer-lock connector.

Packaging: box of 25 individually packaged units.


MD Class IIa Sterilized by ETO

See instructions for use











42, rue Benoit Frachon 93013 - Bobigny Cedex - FRANCE Tel: +33 (0)1 48 10 62 62 Fax: +33(0)1 48 91 22 99 commercial@peters-surgical.com www.peters-surgical.com